EN
CLEANING, DISINFECTION, AND STERILIZATION INSTRUCTIONS FOR STAINLESS STEEL / TITANIUM PRODUCTS
OVERVIEW:

Thorough cleaning and disinfection is indispensable for effective sterilization. The following information refers to a manual procedure with a cleaning/ disinfectant solution prepared by diluting 0.25 fl. oz. of a suitable surgical instrument cleaner in one gallon of de-ionized water. When choosing a cleaning/ disinfectant agent, ascertain that the solution is:

- Suitable for cleaning and disinfecting dental and surgical instruments
- Suitable for ultrasonic cleaning (no foam development)
- Proven efficacy for disinfection
- Compatible with the materials of the products to be cleaned and disinfected
- Free of aldehyde (otherwise blood, secretions, tissue remains, etc. may stick)

CAUTION! FOLLOW THE CLEANING PROCEDURE BELOW TO ACHIEVE EFFECTIVE CLEANING OF DEVICES:

1. Initial disinfection: Immediately after use, immerse all instruments in an industry suitable cleaning/ disinfectant solution for at least 10 minutes.
2. Disassembling: Completely disassemble all instruments consisting of several parts.
3. Preliminary cleaning: Remove coarse contaminants from the instruments. For this purpose use a disposable disinfectant cloth, deionized water, and a soft brush (no metal bristles). Remove residues in areas with difficult access using an appropriate instrument.
4. Combined cleaning and disinfection: Place the instruments for cleaning and disinfection in the cleaning/disinfecting solution for at least 10 minutes, making sure they are completely covered.
 - a. In case of remaining debris, use an ultrasonic unit. Brush off the instruments with a soft brush (no metal bristles) as to completely remove residue.
5. Rinsing: Remove the instruments from the cleaning/ disinfectant solution.
 - a. Next rinse them thoroughly under deionized water.
 - b. Take particular care into areas with limited access, such as channels or grooves.
6. Drying: Dry the instruments. Use a loose, lint-free, disposable disinfectant cloth and compressed dry air.
7. Inspection: Inspect instruments for signs of corrosion, surface damage, chipping, or contamination.
 - a. Remove damaged instruments from use.
 - b. Following inspections, if the instruments are visually contaminated, re-clean and re-rinse the instruments repeating steps 1 through 6 until visual inspection is acceptable.
8. Packaging: Pack the instruments for sterilization promptly using a disposable sterilization package.
 - a. Ensure that the sterilization package is suitable for steam sterilization, and that the products are adequately protected against mechanical damage.

STERILIZATION:

Permissible steam sterilization procedures are fractionated vacuum procedures (with sufficient product drying). Refer to the below validated autoclave parameters for effective sterilization:

1. Sterilization temperature: 270°F (132°C) for 10 minutes
2. The required drying time is 30 minutes (optional)

MANUFACTURER INSTRUCTIONS:

Care for your instruments in accordance with Manufacturer's Instructions.

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