

conFORM™
Resorbable Collagen MEMBRANE

ZERO MEMORY

Memorable Handling.

conFORM™ Resorbable Collagen Membrane, with its **zero memory** characteristics, is an ideal membrane for guiding the healing of bone and surrounding tissue. When hydrated conFORM™ will naturally **conform** to any surgical defect. conFORM's high tensile strength makes it the ideal membrane for suturing, even when hydrated.

conFORM™ is Type 1 bovine collagen and has a predictable resorption time of 3-6 months.



- Fully resorbable
- Easy to handle
- Does not adhere to itself when hydrated
- Zero memory
- 3-6 month predictable resorption

conFORM™

Resorbable Collagen MEMBRANE

conFORM™ RESORBABLE COLLAGEN MEMBRANE is intended for use in oral surgical procedures as a resorbable material for placement to aid in wound healing in the areas of dental implants, bone defects, periodontal defects and ridge augmentation.

ACE conFORM™, with its defined geometry, *in vivo* stability, permeability and mechanical strength, is a white, nonfriable membrane matrix engineered from highly purified Type I Collagen. ACE conFORM™ is packaged sterile, non-pyrogenic, and is for single use. Its resorbability eliminates the need for a second surgery.

FEATURES


- Fully resorbable
- Easy to handle
- Guides healing of bone and surrounding tissue.
- Unique fiber orientation to provide high tensile strength for *in situ* stability of membrane.
- Predictable resorption time (3-6 months).
- Cell occlusive, designed to retard and/or prevent epithelium down growth.

ADVANTAGES

- Does not adhere to itself
- Zero memory
- Sufficient resorption time, 3-6 months, for sustained function in bone repair.
- Higher mechanical strength for membrane stabilization.
- Each piece can be easily trimmed during surgery for a final fit.
- May be placed with either side against the tissue or bone.

TYPICAL USES

- Ridge Augmentation
- Bone Regeneration Around Implants
- Sinus Lift
- Periodontal Defects
- Apical Root Resection
- Extraction Sites



AVAILABLE SIZES

509-1520	15x20mm	\$89.00
509-2030	20x30mm	109.00
509-3040	30x40mm	179.00

All Prices Subject to Change Without Notice