



Once again  
**FASTER** is not  
always **BETTER**

You **CAN'T LOSE**  
with **ACE's RCM<sup>6</sup>**  
**6 Month RESORBABLE**  
**COLLAGEN MEMBRANE**

- **FULLY RESORBABLE**
- **CONTOURS WELL TO DEFECTS**
- **EASILY CUT, SHAPED & APPLIED**
- **COST EFFECTIVE**



**RCM<sup>6</sup>**  
Resorbable Collagen MEMBRANE

To order from  call toll-free **800.441.3100**

[www.acesurgical.com](http://www.acesurgical.com)

ACE Surgical Supply Company, Inc. 1034 Pearl Street, Brockton, MA 02301



Resorbable Collagen MEMBRANE

## RCM<sup>6</sup> RESORBABLE COLLAGEN MEMBRANE

**RESORBABLE COLLAGEN MEMBRANE** is intended for use in oral surgical procedures as a resorbable material for placement to aid in wound healing in the areas of dental implants, bone defects, periodontal defects and ridge augmentation.

**ACE RCM<sup>6</sup>**, with its defined geometry, *in vivo* stability, permeability and mechanical strength, is a white, nonfriable membrane matrix engineered from highly purified Type I Collagen. ACE RCM<sup>6</sup> is packaged sterile, non-pyrogenic, and is for single use. Its resorbability eliminates the need for a second surgery.

### FEATURES

- Guides healing of bone and surrounding tissue.
- Unique fiber orientation to provide high tensile strength for *in situ* stability of membrane.
- Predictable resorption time (26-38 weeks).
- Cell occlusive, designed to retard and/or prevent epithelium down growth.

### ADVANTAGES



- Longer resorption time, 26-38 weeks, for sustained function in bone repair.
- Higher mechanical strength for membrane stabilization.
- Each piece can be easily trimmed during surgery for a final fit.
- May be placed with either side against the tissue or bone.

### TYPICAL USES

- Immediate Implant Loading
- Ridge Augmentation
- Bone Regeneration Around Implants
- Sinus Lift
- Apical Root Resection
- Extraction Sites

### AVAILABLE SIZES

509-6000	15 x 20mm,	.3mm thick	\$77.00
509-3000	20 x 30mm,	.3mm thick	96.00
509-8000	30 x 40mm,	.3mm thick	149.00



### ACE PGA-RAPID ABSORBING\* SUTURES EQUIVALENT TO VICRYL RAPIDE<sup>®</sup>\*\*

Item #	Size	Description	Needle	In.	Qty.	
003-9907	3-0	Undyed Braided PGA-RA	DSM18	18"	12/box	\$76.10
003-9906	4-0	Undyed Braided PGA-RA	DSM18	18"	12/box	76.10
003-9900	4-0	Undyed Braided PGA-RA	HR17 (Taper)	27"	12/box	34.90
003-9904	5-0	Undyed Braided PGA-RA	DSM13	18"	12/box	76.10

\*Maintains 50% tensile strength at 5 days; 0% at 14 days.

\*\*Vicryl<sup>®</sup> and Vicryl Rapide<sup>®</sup> are Registered Trademarks of Ethicon, Inc.

All Prices Subject to Change Without Notice

[www.acesurgical.com](http://www.acesurgical.com)

Fax: 800-583-3150 • Tel: 800-441-3100 • 508-588-3100  
 1034 Pearl Street, P.O. Box 1710, Brockton • MA 02303  
 Email: [info@acesurgical.com](mailto:info@acesurgical.com)