




**conFORM™**  
Resorbable Collagen MEMBRANE

# ZERO MEMORY

## Memorable Handling.

**conFORM™** Resorbable Collagen Membrane, with its zero memory characteristics, is an ideal membrane for guiding the healing of bone and surrounding tissue. When hydrated **conFORM™** will naturally conform to any surgical defect. **conFORM™**'s high tensile strength makes it the ideal membrane for suturing, even when hydrated.

**conFORM™** is Type 1 bovine collagen and has a predictable resorption time of 3-4 months.



- Fully resorbable
- Easy to handle
- Does not adhere to itself when hydrated
- Zero memory
- 3-4 month predictable resorption

# conFORM™

Resorbable Collagen MEMBRANE

**conFORM™ RESORBABLE COLLAGEN MEMBRANE** is intended for use in oral surgical procedures as a resorbable material for placement to aid in wound healing in the areas of dental implants, bone defects, periodontal defects and ridge augmentation.

ACE conFORM™, with its defined geometry, *in vivo* stability, permeability and mechanical strength, is a white, nonfriable membrane matrix engineered from highly purified Type I Collagen. ACE conFORM™ is packaged sterile, non-pyrogenic, and is for single use. Its resorbability eliminates the need for a second surgery.

### FEATURES

- Fully resorbable
- Easy to handle
- Guides healing of bone and surrounding tissue.
- Unique fiber orientation to provide high tensile strength for *in situ* stability of membrane.
- Predictable resorption time (3-4 months).
- Cell occlusive, designed to retard and/or prevent epithelium down growth.

### ADVANTAGES

- Does not adhere to itself
- Zero memory
- Sufficient resorption time, 3-4 months, for sustained function in bone repair.
- Higher mechanical strength for membrane stabilization.
- Each piece can be easily trimmed during surgery for a final fit.
- May be placed with either side against the tissue or bone.

### TYPICAL USES

- Ridge Augmentation
- Bone Regeneration Around Implants
- Sinus Lift
- Periodontal Defects
- Apical Root Resection
- Extraction Sites

### AVAILABLE SIZES

Number	Description	Membrane size	Price
509-1520	conFORM™ Resorbable Collagen Mambrane	15x20mm	\$109.99
509-2030	conFORM™ Resorbable Collagen Mambrane	20x30mm	134.99
509-3040	conFORM™ Resorbable Collagen Mambrane	30x40mm	219.99

