

WHAT YOU AIM FOR IN BONE & MEMBRANE FIXATION

truFIX
A Complete Fixation System

HITS THE MARK EVERY TIME.

truFIX
A Complete Fixation System.

Introducing the *truFIX* System.

The truFIX System from ACE Surgical Supply provides you with everything you need for both bone fixation and membrane stabilization. This system incorporates all the necessary components to pick-up and drive the newly designed ACE truSCREW and truTACK. The truFIX system eliminates the need for multiple systems and unnecessary components, making it user-friendly for your practice.



The newly designed truSCREW, with its aggressive cutting flutes, is the ideal bone screw for the fixation of small bone within the oral and maxillofacial environment. These cutting flutes (in most instances) eliminate the need for any pre-drilling. The patented design of the screw ensures an effortless insertion into all types of bone.


The truTACK makes the stabilization of membranes quick and trouble-free. Our unique tack incorporates a hexagon on its head and threads on its shaft, allowing for easy removal. The truTACK is placed like a tack and removed like a screw, a feature that you will not find in any other system.

The truFIX system is just one more reason that ACE Surgical Supply is your complete dental solution.



www.acesurgical.com

ACE Surgical Supply Company, Inc. • 1034 Pearl Street, Brockton, MA 02301



truFIX™ COMPLETE FIXATION SYSTEM

- You Choose the Screws and Tacks
- truSCREW Developed For Easy Insertion with Maximum Fixation
- True Centered Patented Locking Connection
- True Axial Alignment with Pickup Each Time
- True Self-Drilling Screw
- Patented Removal Sleeve for Disengaging the Screw From the Driver

454-18015 truFIX Complete Fixation System \$899.99

System Includes:

- 454-19001 truFIX Driver Handle
- 454-17001 truTACK Driver Tip
- 454-15201 truSCREW Driver Tip, Contra-Angle
- 454-15001 truSCREW Driver Tip
- 454-17901 truSCREW Driver Removal Sleeve
- 50-27-00045 Pilot Bur, 0.45mm, 27mm Length, Contra-Angle
- 454-40201 1.1mm Pilot Twist Drill
- 454-281001 Step Drill (Lag Screw), .8mm, 1.6mm Steps, Contra-Angle
- 454-50201 truSCREW Packaging Removal Tool
- 20-1577-02 0.88mm External Hex Hand Driver
- 454-18501 truFIX Tray

truSCREW™ 1.2 Thread Diameter

		Head Diameter	Thread Diameter	Thread Length		
454-27202	truSCREW	3HD	1.2TD x 4.5TL	5/pkg	\$109.99	
454-27203	truSCREW	3HD	1.2TD x 6.0TL	5/pkg	109.99	
454-27204	truSCREW	3HD	1.2TD x 7.5TL	5/pkg	109.99	
454-27205	truSCREW	3HD	1.2TD x 9.0TL	5/pkg	109.99	
454-27206	truSCREW	3HD	1.2TD x 10.5TL	5/pkg	109.99	

truSCREW™ 1.5 Thread Diameter

454-27502	truSCREW	3HD	1.5TD x 6.0TL	5/pkg	109.99
454-27503	truSCREW	3HD	1.5TD x 7.5TL	5/pkg	109.99
454-27504	truSCREW	3HD	1.5TD x 9.0TL	5/pkg	109.99
454-27505	truSCREW	3HD	1.5TD x 10.5TL	5/pkg	109.99
454-27506	truSCREW	3HD	1.5TD x 12.0TL	5/pkg	109.99
454-27507	truSCREW	3HD	1.5TD x 13.5TL	5/pkg	109.99
454-27508	truSCREW	3HD	1.5TD x 15.0TL	5/pkg	109.99

truSCREW™ 2.0 Thread Diameter

454-28002	truSCREW	3HD	2.0TD x 6.0TL	5/pkg	109.99
454-28003	truSCREW	3HD	2.0TD x 7.5TL	5/pkg	109.99
454-28004	truSCREW	3HD	2.0TD x 9.0TL	5/pkg	109.99
454-28005	truSCREW	3HD	2.0TD x 10.5TL	5/pkg	109.99
454-28006	truSCREW	3HD	2.0TD x 12.0TL	5/pkg	109.99
454-28007	truSCREW	3HD	2.0TD x 13.5TL	5/pkg	109.99
454-28008	truSCREW	3HD	2.0TD x 15.0TL	5/pkg	109.99

TENTING SCREW 1.5 Thread Diameter

454-25501	Tenting Screw	4HD	1.5TD x 7.0TL	1 each	34.99
454-25502	Tenting Screw	4HD	1.5TD x 11.0TL	1 each	34.99
454-25503	Tenting Screw	4HD	1.5TD x 15.0TL	1 each	34.99

truTACK™

96-003-13	truTACK	2.5HD	3.0TL	10/pkg	115.99
96-003-14	truTACK	2.5HD	5.0TL	10/pkg	115.99

www.acesurgical.com

www.acesurgical.com

Tel: 800-441-3100 • 508-588-3100 • Fax: 800-583-3150
 1034 Pearl Street, Brockton, MA 02301
 Email: info@acesurgical.com

To place an order call toll-free **800.441.3100** 19 | 3-10C